

EDTA with R-Lipoic Acid

Quicksilver Scientific *Etheric Delivery*™ Phospholipid Encapsulation System

4-oz. PET Bottle



Quicksilver's Etheric Delivery™ Phospholipid Encapsulation System brings the power of intravenous therapy into a convenient oral delivery.

The absorption of conventional oral EDTA is only ~5% Nanosphere delivery greatly increases absorption and for some compounds can provide higher *intracellular* delivery than an IV can!

Quicksilver's Etheric Delivery™ Phospholipid Encapsulation System improves upon liposomal technology with smaller, more stable, single-layer spheres made from the highest-grade ingredients available. In addition to exceptional absorption rates, the nanospheres that Quicksilver produces have demonstrated the ability to cross the blood-brain barrier, deposit their

cargo intracellularly and enhance lymphatic circulation of nutrients. Additionally, the phospholipids that compose the liposome shell feed the cell membranes. This ensures the proper function for the absorption of nutrients and the excretion of cellular waste products and toxins.

Quicksilver Scientific's Etheric Delivery™ EDTA w/ Sodium R-Lipoic Acid is the most absorbable form of oral EDTA. EDTA has broad metals detoxification ability, being especially for lead and R-Lipoic acid (as sodium R-Lipoate) has an exceptionally well-documented ability to upregulate the glutathione system via the Nrf2 nuclear transcription pathway. Combined this EDTA/R-Lipoate with Etheric Delivery Glutathione and IMD, and you have an incredible, completely-oral broad-spectrum metal detoxification system that not only eliminates the dependence on IV's, but repairs the underlying problem in the detoxification system that led to metals accumulation in the first place. Increased lead excretion is comparable to IV and continues for several days after a single bolus dose.

Production Technology Features:

- Tightly controlled small sizing, verified by Laser Dynamic Light Scattering on each batch.
- Mean size 50-100nm
- Flavor suitable for intra-buccal delivery, or easy water-dispersability when desired

General Suggested Usage: For general metal detoxification, take 1 tsp (210 mg EDTA and 33 mg of R-Lipoate) of product daily either straight by mouth or diluted into water. Alternately, it can be taken as a 5-tsp dose (1050mg EDTA and 165mg R-Lipoate)

1-2X per week. Be sure to use adequate mineral supplementation, especially zinc, when using this product.

SUPPLEMENT FACTS		
Serving Size: 1 tsp (5 mL)	Amount	% Daily
Servings per container: 24	Per Serving	Value
Sodium	68 mg	2.8
EDTA (Disodium EDTA)	210 mg	*
R-Lipoic Acid	33 mg	*
Phosphatidyl Choline (Soy Lecithin - 65% linoleic acid)	450 mg	*

* Daily Value not established

Other Ingredients: Water, Glycerin, Ethanol, sodium hydroxide

